

THE CUTTING EDGE

KEIRSURGICAL

October/November 2009

Keir Surgical's Perforated Coagulating Resector

Keir Surgical's Perforated Coagulating Resectors are a unique and highly regarded product for surgical intervention of Benign Prostatic Hyperplasia (BPH). Trans-Urethral Resection of the Prostate (TURP) offers the Urologist the most effective approach to removal of prostate tissue. Keir Surgical's Perforated Coagulating Resector gives the surgeon the perfect balance of cut and coagulation. Cutting is optimized because the same thickness as a standard cutting loop (.012 inches) but the tip is wider and perforated to enhance coagulation.



Some additional advantages of the Keir Surgical Perforated Coagulating Resectors are:

- *The best cutting outside of a traditional cutting loop*
- *Maintains the thin "cutting loop profile" to allow surgical debulking of the prostate*
- *Perforated loop creates "current density", enhancing coagulation by increasing the depth of energy penetration*
- *Broad loop sustains contact with the tissue longer, allowing deeper penetration of energy and increasing coagulation*

Available in Karl Storz, ACMI, and Olympus-compatible models.

LUKENS® Bone Wax



Lukens® Bone Wax has a new formulation! It is softer, easier to handle, form and apply than ever before. The peel-open foil packet (packed 12 to a box) is specifically designed for easy dispensing and provides optimum safety and convenience for use in the operating room.

Bone Wax is used to achieve local hemostasis of bone by acting as a mechanical sealant (tamponade). It does not act biochemically and is minimally resorbable.

Lukens™ Bone Wax is used to control bleeding on bone surface in:

- Orthopedic Surgery
- Thoracic Surgery
- Neurosurgery
- Oral & Maxillofacial Surgery

Call Keir Surgical today for a free sample!

KS

KEIRSURGICAL

126-408 East Kent Avenue South
Vancouver, BC V5X 2X7

T: 800.663.4525 | F: 604.261.9549

info@keirsurgical.com | www.keirsurgical.com

Please contact your Keir Surgical Sales Representative for more information.

THE CUTTING EDGE


KEIRSURGICAL

October/November 2009

ACCESS Surgical Bariatric Instruments

Access Surgical Bariatric instruments are reusable and offer the precision, performance and durability that top surgeons demand. These instruments are specifically designed for the Gastric Bypass Procedure (circular anvil grasper and perforator) and the Jones Gastric Banding Procedure (hook, pusher, and band-passer).

With unique handles that fit comfortably in the surgeon's palm, these instruments provide optimal comfort and ergonomic ease. Access Surgical Bariatric Instruments are designed with a flush port that speeds effective cleaning of inner lumens. These instruments also carry a lifetime warranty against defects on material and workmanship.



Did You Know?

- Keir Surgical has Canadian exclusivity with more than 15 suppliers including Symmetry Medical (Riley Medical, Access Surgical, and SSI Ultra), Rudolf Medical, BFW Inc., and Mediflex.
- Keir Surgical is active in the communities in which we live through frequent donations and volunteer actions; we are proud sponsors of the Vancouver Symphony Orchestra and have donated over \$1500 to Canadian charities through employee donation matching programs in the past 15 months.
- Keir Surgical has six instructional videos on its website, addressing commonly asked product questions. To view these videos, please visit www.keirsurgical.com/education_videos.html.
- Upcoming conferences:

CAOMS (Canadian Oral Maxillofacial) – October 15-17 in Toronto, ON

ORNAA (Alberta OR Nurses) – October 21-24 in Red Deer, AB

ORNAHD (Hamilton and District OR Nurses) – October 24 in Niagara on the Lake, ON

AGS conference (Quebec Sterilization Managers) – October 29 & 30 in Québec City, QC

OAGS (Ontario General Surgeons) – November 7 in Toronto, ON

How are we doing?

Take the new Keir Surgical Customer Satisfaction Survey.

www.keirsurgical.com/satisfactionsurvey

KS

KEIRSURGICAL

126-408 East Kent Avenue South
Vancouver, BC V5X 2X7

T: 800.663.4525 | F: 604.261.9549

info@keirsurgical.com | www.keirsurgical.com

Please contact your Keir Surgical Sales Representative for more information.